G MIXER

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GANDO PRODUCTS

G MIXER

VACUUM MIXER FOR DENTAL LABS; THE NEW GENERATION OF PLASTER MIXING UNDER VACUUM

IN MODERN DENTAL TECHNOLOGY, THE QUALITY OF MATERIALS AND THE PRECISION OF THEIR PREPARATION HAVE A DIRECT IMPACT ON THE FINAL OUTCOME OF DENTAL PROSTHESES, IMPRESSIONS, AND LABORATORY MODELS. ONE OF THE MOST ESSENTIAL TOOLS IN THIS PROCESS IS THE DENTAL VACUUM MIXER, AN ADVANCED DEVICE SPECIFICALLY DESIGNED FOR MIXING DENTAL PLASTERS UNDER VACUUM CONDITIONS, EFFECTIVELY PREVENTING THE FORMATION OF AIR BUBBLES IN THE MIXTURE.

THE PRESENCE OF AIR BUBBLES IN DENTAL PLASTER HAS ALWAYS BEEN A MAJOR CHALLENGE IN DENTAL LABORATORIES. EVEN THE SMALLEST BUBBLES CAN REDUCE THE STRENGTH OF THE PLASTER, CREATE SURFACE DEFECTS IN MODELS, DECREASE IMPRESSION ACCURACY, AND ULTIMATELY INCREASE ERRORS IN PROSTHESIS FABRICATION. BY CREATING A VACUUM ENVIRONMENT DURING MIXING, THE VACUUM MIXER REMOVES TRAPPED AIR FROM THE MIXTURE, RESULTING IN A HOMOGENEOUS, DENSE, AND BUBBLE-FREE COMPOUND.

DIFFERENCES BETWEEN THE GANDO VACUUM MIXER AND CONVENTIONAL DEVICES

THE NEW-GENERATION VACUUM MIXER, DESIGNED BY ADVANCED ENGINEERING TEAMS, HAS SIGNIFICANT IMPROVEMENTS COMPARED TO OLDER MODELS AND EXISTING MARKET DEVICES. SOME OF THESE DIFFERENCES INCLUDE:

HIGHER VACUUM POWER: THE VACUUM SYSTEM IN THIS DEVICE IS ENGINEERED TO MINIMIZE TRAPPED AIR IN THE PLASTER TO THE LOWEST POSSIBLE LEVEL, WHEREAS MANY OLDER DEVICES ARE UNABLE TO MAINTAIN A STABLE VACUUM.

ERGONOMIC AND USER-FRIENDLY DESIGN: UNLIKE TRADITIONAL MODELS THAT WERE OFTEN COMPLICATED TO USE, THIS DEVICE FEATURES A DIGITAL CONTROL PANEL, A CLEAR DISPLAY, AND PRE-SET OPERATING PROGRAMS. IT ALSO INCLUDES THREE WORKING ALGORITHMS: PER VACUUM, PER MIX, AND POST VACUUM.

LOW-NOISE AND LOW-VIBRATION MOTOR: ONE OF THE COMMON ISSUES IN OLDER DEVICES WAS EXCESSIVE NOISE AND VIBRATION DURING OPERATION. THIS NEW VACUUM MIXER, EQUIPPED WITH A HIGH-PERFORMANCE MOTOR (BÜHLER) AND A DYNAMIC BALANCING SYSTEM, ENSURES SMOOTH AND STABLE OPERATION.

HIGH DURABILITY AND LONG LIFESPAN: THE USE OF HIGH-QUALITY MATERIALS, CORROSION-RESISTANT ALLOYS, AND INDUSTRIAL-GRADE COMPONENTS ENSURES LONG-TERM PERFORMANCE WITHOUT QUALITY DEGRADATION.

EASY CLEANING AND MAINTENANCE: THE MODULAR DESIGN OF THE MIXING CHAMBER ALLOWS FOR QUICK CLEANING AND DISINFECTION, WHICH IS ESPECIALLY IMPORTANT IN DENTAL LABORATORY ENVIRONMENTS.

GANDO PRODUCTS

G-MIXER GANDO; A PROGRAMMABLE NEXT-GENERATION VACUUM MIXING SYSTEM

THE G-MIXER GANDO IS AN ADVANCED PROGRAMMABLE DENTAL DEVICE WITH THE ABILITY TO STORE CUSTOMIZED OPERATING PROGRAMS IN ITS INTERNAL MEMORY. USERS CAN SAVE UP TO TEN DIFFERENT PROGRAMS, ALLOWING FOR PRECISE CONTROL OVER VARIOUS MIXING .PROTOCOLS BASED ON DIFFERENT MATERIALS AND CLINICAL REQUIREMENTS
IN ADDITION, THE DEVICE ALSO SUPPORTS SAVING KEY PARAMETERS SUCH AS MIXING TIME AND MOTOR SPEED, ENABLING CONSISTENT AND REPEATABLE PERFORMANCE FOR .LABORATORY WORKFLOWS

WHAT TRULY DIFFERENTIATES THIS VACUUM MIXER FROM CONVENTIONAL SYSTEMS ARE ITS :THREE INNOVATIVE FUNCTIONAL ELEMENTS

PER VACUUM

,THIS FUNCTION ACTIVATES THE VACUUM SYSTEM BEFORE THE MIXING PROCESS BEGINS INCREASING THE VACUUM LEVEL TO ITS MAXIMUM WITHIN A DEFINED TIME RANGE (FROM 1 ,TO 9 SECONDS). THIS ENSURES THAT THE CHAMBER IS FULLY EVACUATED PRIOR TO MIXING .MINIMIZING INITIAL AIR ENTRAPMENT

PER MIX

THIS MODE SIMULATES THE ACTION OF MANUAL SPATULATION. WITH THIS FEATURE, THERE IS NO NEED TO MANUALLY MIX PLASTER AND LIQUID WITH A SPATULA BEFORE PLACING THE CONTAINER INTO THE DEVICE. THE PER MIX FUNCTION PERFORMS THIS INITIAL BLENDING .AUTOMATICALLY AND CAN BE ADJUSTED BETWEEN 0 TO 10 SECONDS

POST VACUUM

THE THIRD FUNCTION IS ACTIVATED AFTER THE MIXING PROCESS IS COMPLETED. BY ,APPLYING A STRONG FINAL VACUUM PHASE, IT REMOVES ANY REMAINING MICRO-BUBBLES .ENSURING A DENSE, UNIFORM, AND COMPLETELY AIR-FREE MIXTURE

:IF YOU WANT, I CAN NEXT

TURN THIS INTO A PROFESSIONAL BROCHURE LAYOUT

OR WRITE A HIGH-IMPACT SALES SCRIPT FOR VIDEO NARRATION

OR SHORTEN IT INTO A WEBSITE LANDING PAGE VERSION

THE G-MIXER IS EQUIPPED WITH TWO TRANSPARENT MIXING CUPS IN 200 ML AND 500 ML SIZES, ALLOWING FLEXIBILITY FOR DIFFERENT WORKING VOLUMES AND MATERIAL REQUIREMENTS IN DENTAL LABORATORY PROCEDURES.

THE UPPER MIXING SECTION IS MADE OF HIGH-QUALITY ALUMINUM, PROVIDING EXCELLENT DURABILITY, A LIGHTWEIGHT STRUCTURE, AND A POLISHED, PREMIUM APPEARANCE. THIS DESIGN ENSURES LONG-TERM PERFORMANCE WHILE MAINTAINING A PROFESSIONAL FINISH.

THE SHAFT AND MIXING COMPONENTS ARE MANUFACTURED FROM 304 STAINLESS STEEL, WHICH OFFERS HIGH RESISTANCE TO CORROSION AND WEAR. THIS ENSURES LONG-LASTING STABILITY, HYGIENE, AND RELIABILITY EVEN UNDER CONTINUOUS USE IN DEMANDING LABORATORY ENVIRONMENTS.

APPLICATIONS OF THE GANDO VACUUM MIXER IN DENTAL LABORATORY WORK

THE GANDO VACUUM MIXER IS WIDELY USED ACROSS A BROAD RANGE OF DENTAL LABORATORY PROCEDURES, INCLUDING:

- PREPARATION OF GYPSUM MODELS FOR FIXED AND REMOVABLE PROSTHESES
- INVESTING (POURING) CYLINDERS FOR INJECTION SYSTEMS USED IN LAMINATE VENEERS
- PRECISE IMPRESSION AND MODEL PREPARATION FOR IMPLANT CASES
- FABRICATION OF WORKING MODELS FOR CROWNS, VENEERS, AND BRIDGES

BY USING A VACUUM MIXER, THE DENTAL TECHNICIAN CAN ENSURE THAT THE FINAL MODEL IS COMPLETELY FREE OF AIR BUBBLES AND STRUCTURAL INCONSISTENCIES. THIS RESULTS IN HIGHER ACCURACY, IMPROVED MATERIAL STRENGTH, AND A SIGNIFICANTLY MORE PROFESSIONAL AND RELIABLE FINAL OUTCOME.

GANDO PRODUCTS

AUTOMATIC START FUNCTION

THE G-MIXER IS EQUIPPED WITH TWO MICRO-SWITCHES LOCATED AT ITS BASE, ENABLING A FULLY AUTOMATIC OPERATION MODE. BY SIMPLY PLACING THE MIXING CUP UNDER THE MIXING UNIT, THE DEVICE AUTOMATICALLY DETECTS IT AND STARTS THE PROCESS WITHOUT THE NEED TO PRESS ANY BUTTONS

THIS FEATURE PROVIDES A MORE EFFICIENT, FASTER, AND USER-FRIENDLY WORKFLOW, REDUCING MANUAL INTERACTION AND IMPROVING CONSISTENCY IN LABORATORY OPERATIONS

FLEXIBLE INSTALLATION DESIGN

THE G-MIXER OFFERS A FLEXIBLE INSTALLATION SYSTEM THAT ALLOWS USERS TO OPTIMIZE THEIR WORKSPACE ACCORDING TO THEIR LABORATORY NEEDS. BY REMOVING THE BASE STAND, THE DEVICE CAN BE MOUNTED ON THE WALL, SIGNIFICANTLY REDUCING OCCUPIED BENCH SPACE AND IMPROVING WORKFLOW EFFICIENCY

ALTERNATIVELY, THE UNIT CAN ALSO BE USED IN A STANDARD CONFIGURATION BY INSTALLING IT ON ITS DEDICATED BASE. THIS DUAL INSTALLATION OPTION PROVIDES GREATER ADAPTABILITY FOR DIFFERENT LABORATORY ENVIRONMENTS AND SPACE LIMITATIONS

GANDO PRODUCTS

WHY INVESTING IN A VACUUM MIXER IS WORTH IT

ANY DENTAL LABORATORY AIMING TO IMPROVE WORK QUALITY, REDUCE ERRORS, SAVE TIME, AND MINIMIZE THE COST OF REWORK SHOULD CONSIDER USING A STANDARD VACUUM MIXER. THIS DEVICE NOT ONLY SIGNIFICANTLY ENHANCES THE QUALITY OF THE FINAL OUTPUT, BUT ALSO MAKES THE ENTIRE WORKFLOW MORE PROFESSIONAL AND EFFICIENT, ULTIMATELY INCREASING SATISFACTION FOR BOTH DENTISTS AND PATIENTS.

IF YOU ARE LOOKING TO UPGRADE THE QUALITY OF YOUR LABORATORY WORK AND REDUCE PROCEDURAL ERRORS, INVESTING IN THE GANDO VACUUM MIXER IS A SMART AND FORWARD-THINKING DECISION THAT SUPPORTS LONG-TERM GROWTH AND PRECISION IN DENTAL FABRICATION.

WEIGHT	12 kg
GLASS	200 ml & 500 ml
VOLTAGE	220 AC
POWER	160W

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